

# Compounding Formulation Record for: FIRST®-Lansoprazole 3 mg/mL- Suspension Compounding Kit

## Non-FDA Approved Compounding Kit: For Pharmacist Use Only

Patient Name: \_\_\_\_\_

Date Prepared: \_\_\_\_\_

Rx: \_\_\_\_\_ - \_\_\_\_\_

Batch #: \_\_\_\_\_ Beyond Use Date: 30 Days Good For Use Until: \_\_\_\_/\_\_\_\_/\_\_\_\_

**CLASS: RX**      **DOSAGE FORM:** Suspension      **ROUTE OF ADMINISTRATION:** Oral

Ingredient Name	Dispensed Quantity	Product Code	Manufacturer	Lot #	Exp. Date of Ingredient	Prepared By:	Verified By:
Bulk Lansoprazole USP Powder		65628-080-	Azurity				
FIRST® - Suspension		65628-080-	Azurity				

Prior to compounding, store FIRST®-Lansoprazole Compounding Kit at refrigerated temperature 2°C - 8°C (36°F - 46°F)  
Store the final compounded formulation at refrigerated temperature 2°C - 8°C (36°F - 46°F)

FIRST®-Lansoprazole 3 FL OZ (90 mL) Kit (65628-080-03) contains 0.27 g Bulk Lansoprazole USP Powder + 90 mL FIRST® Suspension

FIRST®-Lansoprazole 5 FL OZ (150 mL) Kit (65628-080-05) contains 0.45 g Bulk Lansoprazole USP Powder + 150 mL FIRST® Suspension

FIRST®-Lansoprazole 10 FL OZ (300 mL) Kit (65628-080-10) contains 0.9 g Bulk Lansoprazole USP Powder + 300 mL FIRST® Suspension

Total Quantity	mL						
----------------	----	--	--	--	--	--	--

### SEE PHARMACIST COMPOUNDING INSTRUCTIONS FOR ADDITIONAL INFORMATION.

**COMPOUNDING DIRECTIONS:**

1. FIRST®- Lansoprazole Compounding Kit contains pre-measured lansoprazole powder and FIRST® - Suspension.
2. Hold the neck of the bottle containing the lansoprazole powder and tap the bottom edges on a hard surface to loosen the powder. Remove the cap from the bottle.
3. Tap the top of the induction seal liner to loosen any powder which may have adhered to the liner. Carefully and slowly peel back the inner foil seal liner from the bottle. Using the enclosed tool, scrape any lansoprazole powder from the seal into the bottle.
4. Shake the FIRST® - Suspension bottle for a few seconds. Open the suspension bottle and empty about one-third the contents of the FIRST®- Suspension bottle **into the lansoprazole powder bottle** Replace the cap and shake the lansoprazole powder bottle for approximately 30 seconds. Shake vigorously.
5. Empty another one-third of the contents of FIRST®-Suspension **into the lansoprazole powder bottle**. Replace the cap and shake the lansoprazole powder bottle for approximately 30 seconds. Shake vigorously.
6. Empty the remaining FIRST®- Suspension **into the lansoprazole powder bottle**. **Allow suspension to drain for 10 seconds**. Replace the cap and shake the lansoprazole powder bottle for approximately 30 seconds. Shake vigorously.
7. **Shake the bottle well before each use.**
8. If appropriate, dispense enclosed adapter cap and oral syringe with accompanying instructions for use.

**WARNING: ADAPTER CAP IS NOT CHILD-RESISTANT**

**RECOMMENDED AUXILIARY LABELS:**

- For oral use only
- Protect from light
- Avoid contact with eyes
- Protect from freezing
- Keep container tightly closed
- The beyond-use date of the compounded product, as dispensed, when stored at refrigerated temperature is **not later than 30 days**.
- Keep out of the reach of children

**REFERENCE:**

FIRST® - LANSOPRAZOLE 3 MG/ML [PHARMACIST COMPOUNDING INSTRUCTIONS]  
Woburn, MA: Azurity Pharmaceuticals, Inc. 2024  
**PART CMP-688 | REV. 03**

